

Code No.01-2337

SPECIFICATION

Amidosulfuric Acid

TraceSure(R)

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay (HOSO ₂ NH ₂)	min.99.90%(mass/mass)
Water-insoluble matter	max.0.01%(mass/mass)
Loss on drying at 50°C	max.0.1%(mass/mass)
Residue after ignition	max.0.005%(mass/mass)
Chloride (Cl)	max.5mg/kg
Nitrate (NO ₃)	max.5mg/kg
Copper (Cu)	max.1mg/kg
Lead (Pb)	max.1mg/kg
Iron (Fe)	max.2mg/kg

JIS K 8005:2016

*:Additional test performed by Wako but not required by JIS