

Code No.01-2370

## SPECIFICATION

### Agarose LE (Low electroendosmosis) for Electrophoresis

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown, powder
Solubility in water	to pass test
Gel point	to pass test
Electroendosmosis (-Mr)	0.10~0.15
Loss on drying at 105°C	max.10.0%
Sulfate (SO <sub>4</sub> )	max.0.4%
Jelly strength (1.5%)	min.1200g/cm <sup>2</sup>