

SPECIFICATION

Albumin, from Bovine Serum, Fraction V, pH7.0 for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White ~ brown or brownish green, crystals ~
Solubility in water	crystalline powder to pass test
Loss on drying over silica gel in vacuum	max.5.0%
pH of a 100g/L solution at 25°C	6.5~7.5
Heavy metals (as Pb)	max.0.005%
Bacterial endotoxins	less than 3EU/mg
Electrophoresis test (SDS-PAGE)	to pass test