

Code No.01-2765

SPECIFICATION

Albumin, Human, recombinant expressed in yeast  
for Cell Culture

REQUIREMENT	SPECIFICATION
Appearance	White ~ yellowish brown, crystals ~ powder or mass
Solubility in water	to pass test
Loss on drying over silica gel in vacuum	under investigation
pH of a 20g/L solution at 25°C	under investigation
Bacterial endotoxins	less than 5EU/g
Assay (HPLC)	min.98.0%

Source:Yeast expressed human serum albumin