

Code No.01-2838

SPECIFICATION

Albumin, from Bovine Serum, Fatty Acid/Protease Free, Low IgG for Biochemistry

| REQUIREMENT | SPECIFICATION |
|---------------------------------|---|
| Appearance | Slightly greenish yellow ~ brown, crystals ~ crystalline powder to pass test max.5.0% 6.5~7.5 under investigation to pass test to pass test min.98% |
| Solubility in water | |
| Loss on drying at 135°C | |
| pH of a 100g/L solution at 25°C | |
| IgG | |
| Fatty acid | |
| Protease | |
| Assay (electrophoresis) | |