

Code No.01-2838

SPECIFICATION

Albumin, from Bovine Serum, Fatty Acid/Protease Free, Low IgG for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	Slightly greenish yellow ~ brown, crystals ~
Solubility in water	crystalline powder
Loss on drying at 135°C	to pass test
pH of a 100g/L solution at 25°C	max.5.0%
IgG	6.5~7.5
Fatty acid	under investigation
Protease	to pass test
Assay (electrophoresis)	to pass test
	min.98%