

Code No.02-0342

SPECIFICATION

 (\pm) -1-Methylpropylamine

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ nearly colorless, clear liquid
Solubility in water	to pass test
Density (20°C)	0.720~0.726g/ml
Refractive index 20°C	1.391~1.395
Specific rotation $[\alpha]_{D20}(\text{neat})$	-0.1~+0.1°
Water	max.0.5%
Assay (capillary GC)	min.99.0%