

Code No.02-0550

SPECIFICATION

Berberine Chloride n-Hydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Yellow, powder
Solubility in water	to pass test
Water	max.16.5%
Residue after ignition (as sulfate)	max.0.1%
Acidity (as HCl)	max.0.4%
Sulfate (SO ₄)	max.0.05%
Heavy metals (as Pb)	max.0.003%
Assay (calculated on the anhydrous basis)	min.97.0%