

Code No.02-0827

## SPECIFICATION

## Bis(2-butoxyethyl) Ether

Wako Special Grade

| REQUIREMENT                       | SPECIFICATION                              |
|-----------------------------------|--|
| Appearance                        | Colorless ~ nearly colorless, clear liquid |
| Density (20°C)                    | 0.881~0.887g/ml                            |
| Refractive index 20°C             | 1.422~1.425                                |
| Water                             | max.0.05%                                  |
| Acidity (as CH <sub>3</sub> COOH) | max.0.02%                                  |
| Assay (capillary GC)              | min.98.0%                                  |