

Code No.02-1076

## SPECIFICATION

### Bismuth(III) Subgallate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Yellow powder
Solubility in sodium hydroxide solution	to pass test
Loss on drying at 105°C	max.6.0%
Chloride (Cl)	max.0.05%
Copper (Cu)	max.0.02%
Lead (Pb)	max.0.02%
Iron (Fe)	max.0.03%
Assay (as Bi) (after drying)	47.0~51.0%