

Code No.02-1161

## SPECIFICATION

### Benzyl Methacrylate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear liquid
Solubility in ethanol	under investigation
Density (20°C)	1.036~1.044g/ml
Refractive index 20°C	1.511~1.516
Water	max.0.1%
Acidity [as CH <sub>2</sub> :C(CH <sub>3</sub> )COOH]	max.0.1%
Assay (capillary GC)	min.98.0%