

Code No.02-1836

SPECIFICATION

t-Butyl Methyl Ether, Super Dehydrated
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.739~0.744g/mL
Water	max.0.001%
Residue after evaporation	max.10ppm
Acidity (as CH ₃ COOH)	max.0.001%
Peroxides (as H ₂ O ₂)	max.5ppm
Assay (capillary GC)	min.99.5%