

Code No.03-0044

SPECIFICATION

Calcium Chloride Dihydrate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystalline powder
Solubility in water	to pass test
pH of a 50g/l solution at 25 °C	4.5~9.5
Nitrate (NO ₃)	max.0.003%
Sulfate (SO ₄)	max.0.01%
Heavy metals (as Pb)	max.0.001%
Barium (Ba)	max.0.006%
Arsenic (As)	max.3 ppm
Iron (Fe)	max.0.001%
Ammonium (NH ₄)	max.0.005%
Assay	93.0~103.0%