

Code No.03-0124

SPECIFICATION

Carbon Disulfide

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (CS ₂) (GC)	min.99.0%(mass/mass)
Water	max.0.02%(mass/mass)
Residue after evaporation	max.0.002%(mass/mass)
Acidity	to pass test
Sulfate (SO ₄)	max.3ppm(mass/mass)
Iodine-consuming substances	to pass test

JIS K 8732:2011

(*:Additional test performed by Wako but not required by JIS)