

SPECIFICATION

Carbon Tetrachloride

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	1.592~1.598g/ml
Refractive index 20°C	1.459~1.461
Water	max.0.02%
Residue after evaporation	max.0.001%
Acidity (as HCl)	max.0.001%
Chloride (Cl)	max.1ppm
Free chlorine	to pass test
Carbon disulfide	to pass test
Iodine-consuming substances	to pass test
Substances darkened by sulfuric acid	to pass test
Suitability for use with dithizone	to pass test
Assay (capillary GC)	min.99.5%