

Code No.03-0300

SPECIFICATION

Cholesterol

Wako Special Grade

| REQUIREMENT | SPECIFICATION |
|--|--|
| Appearance | White ~ slightly yellow, crystals ~ powder |
| Solubility in ethanol | to pass test |
| Melting point | 146~149°C |
| Specific rotation [α] _{D25} (c=2, 1,4-dioxane) (after drying) | -34~-39° |
| Loss on drying at 105°C | max.0.2% |
| Acidity | to pass test |
| Assay (capillary GC) | min.90.0% |