

SPECIFICATION

Cobalt(II) Acetate Tetrahydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Reddish purple, crystals ~ crystalline powder
Solubility in dil.acetic acid	to pass test
pH of a 50g/L solution at 25°C	min.5.0
Chloride (Cl)	max.0.003%
Nitrate	to pass test
Sulfate (SO ₄)	max.0.005%
Sodium (Na)	max.0.02%
Potassium (K)	max.0.02%
Copper (Cu)	max.0.005%
Calcium (Ca)	max.0.02%
Zinc (Zn)	max.0.03%
Lead (Pb)	max.0.005%
Iron (Fe)	max.0.005%
Nickel (Ni)	max.0.3%
Assay	min.99.0%