

Code No.03-0424

## SPECIFICATION

### Copper(II) Nitrate Trihydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Solubility in water	to pass test
Solubility in ethanol	to pass test
Chloride (Cl)	max.0.002%
Sulfate (SO <sub>4</sub> )	max.0.005%
Sodium (Na)	max.0.005%
Potassium (K)	max.0.005%
Calcium (Ca)	max.0.005%
Zinc (Zn)	max.0.005%
Iron (Fe)	max.0.003%
Nickel (Ni)	max.0.002%
Assay	99.0~104.0%