

Code No.03-0854

## SPECIFICATION

## p-Isopropylbenzaldehyde

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.977~0.983g/ml
Refractive index 20°C	1.527~1.533
Water	max.0.2%
Acidity [as (CH <sub>3</sub> ) <sub>2</sub> CHC <sub>6</sub> H <sub>4</sub> COOH]	max.0.5%
Assay (capillary GC)	min.95.0%