

SPECIFICATION

Carbon Tetrachloride for HPLC

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Refractive index 20°C	1.459~1.461
Absorbance (265nm)	max.0.70
Absorbance (270nm)	max.0.35
Absorbance (280nm)	max.0.06
Absorbance (290nm)	max.0.02
Absorbance (295~400nm)	max.0.01
Water	max.0.01%
Acidity (as HCl)	max.2ppm
Free chlorine	to pass test
Fluorescence test	to pass test
Assay (capillary GC)	min.99.5%