

Code No.03-0921

SPECIFICATION

Cyclohexyl Methacrylate

Wako Special Grade

| REQUIREMENT | SPECIFICATION |
|---|---|
| Appearance | Colorless ~ slightly yellow, clear liquid |
| Density (20°C) | 0.963~0.970g/mL |
| Refractive index 20°C | 1.456~1.461 |
| Water | max.0.5% |
| Acidity [as CH ₂ :C(CH ₃)COOH] | max.0.1% |
| Assay (capillary GC) | min.97.0% |