

Code No.03-1074

## SPECIFICATION

## DL-p-Chlorophenylalanine

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown, crystalline powder ~
Solubility in hydrochloric acid	powder to pass test
Specific rotation [ $\alpha$ ]D20(c=0.5,dil.HCl)	under investigation
Residue after ignition (as sulfate)	max.0.1%
Thin layer chromatography test (1)	to pass test
Thin layer chromatography test (2)	to pass test
Assay	min.99.0%