

Code No.03-1313

## SPECIFICATION

Choline Chloride  
for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Absorbance (500g/L, 230nm)	max.0.10
Absorbance (500g/L, 260nm)	max.0.20
Loss on drying at 110°C	max.1.0%
pH of a 50g/L solution at 25°C	4.0~7.0
Heavy metals (as Pb)	max.0.001%
Ammonium (NH <sub>4</sub> )	max.0.01%
Assay	min.98.0%