

Code No.03-1953

SPECIFICATION

Chlormadinone Acetate Standard for High Performance Liquid Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale yellow, crystals ~ powder
Solubility in acetonitrile	to pass test
Specific rotation [α]D ₂₀ (c=2,CH ₃ CN) (after drying)	-10~-14°
Loss on drying over phosphorus(V) oxide in vacuum	max.0.5%
Foreign steroids	to pass test
Assay (HPLC)	min.99.0%