

Code No.03-1969

SPECIFICATION

Cephalonium Dihydrate Standard

for High Performance Liquid Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale yellowish brown,
Solubility in dimethyl sulfoxide	crystalline powder ~ powder
Specific rotation [ $\alpha$ ]D20(c=0.5,DMSO,30min)	to pass test
(calculated on the dehydrous basis)	
Water	under investigation
Assay (HPLC)	6.5~8.0%
	min.96.0%