

SPECIFICATION

Carbazochrome Sodium Sulfonate Trihydrate
for Crude Drugs Determination

REQUIREMENT	SPECIFICATION
Description (JP method)	Reddish yellow, crystals or crystalline powder to pass test
Identification test (JP method)	to pass test
pH (JP method)	5.0~6.0
Purity (JP method)	-
(1) Clarity and color of solution (JP method)	to pass test
(2) Heavy metals (JP method)	max.20ppm
(3) Analogous substances (JP method)	to pass test
Water	14.0~15.0%
Assay (calculated on the dehydrous basis)	min.99.0%