

SPECIFICATION

Calcium Chloride Dihydrate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 50g/L solution at 25°C	4.5~9.5
Nitrate (NO ₃)	max.0.003%
Sulfate (SO ₄)	max.0.01%
Copper (Cu)	max.2ppm
Barium (Ba)	max.0.006%
Lead (Pb)	max.0.001%
Arsenic (As)	max.3ppm
Iron (Fe)	max.0.001%
Ammonium (NH ₄)	max.0.005%
Assay	93.0~103.0%