

SPECIFICATION

CultureSure™ 3,3',5-Triiodo-L-thyronine Sodium Salt (T3)
for Cell Culture

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale yellowish brown or reddish brown, crystalline powder ~ powder
Solubility in dimethyl sulfoxide	to pass test
Specific rotation [α] _{D25}	+18.0~+22.0°
Water	measured value
Bacterial endotoxins	less than 0.25EU/mg
Mycoplasma test	to pass test
Assay (HPLC)	min.98.0%