

Code No.04-0276

SPECIFICATION

N,N-Dimethylaniline

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear liquid
Density (20°C)	0.954~0.958g/mL
Refractive index 20°C	1.556~1.560
Water	max.0.2%
Aniline (C ₆ H ₅ NH ₂)	max.0.03%
N-Methylaniline (C ₆ H ₅ NHCH ₃)	max.0.3%
Assay (capillary GC)	min.99.0%