

Code No.04-0277

SPECIFICATION

N,N-Dimethylaniline

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear liquid
Density (20°C)	0.953~0.959g/mL
Refractive index 20°C	1.556~1.560
Water	max.0.5%
Aniline (C ₆ H ₅ NH ₂)	max.0.03%
N-Methylaniline (C ₆ H ₅ NHCH ₃)	max.0.5%
Assay (capillary GC)	min.99.0%