

## SPECIFICATION

## N,N-Dimethylformamide

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Description	Colorless clear liquid
Assay [HCON(CH <sub>3</sub> ) <sub>2</sub> ] (GC)	min.99.5%(mass/mass)
Appearance	max.10 Hazen units
Solubility in water	to pass test
Solubility in ethanol	to pass test
Density (20°C)	0.947~0.952g/ml
Refractive index 20°C	1.428~1.432
Water	max.0.1%(mass/mass)
Residue after evaporation	max.0.005%(mass/mass)
Acidity (as HCOOH)	max.0.003%(mass/mass)
Alkalinity (as NH <sub>3</sub> )	max.0.002%(mass/mass)
Iron (Fe)	max.5ppm(mass/mass)

JIS K 8500:2007(\*:Additional test performed by Wako but not required by JIS)