

Code No.04-0424

SPECIFICATION

2,2'-Bipyridyl

Guaranteed Reagent

| REQUIREMENT | SPECIFICATION |
|--|--|
| * Appearance | White ~ slightly brown, crystals ~ crystalline powder |
| Assay (C ₁₀ H ₈ N ₂) | min.99.5%(mass/mass) |
| Solubility in ethanol | to pass test |
| Melting point | 69.5~72.0°C |
| Residue after ignition (as sulfate) | max.0.05%(mass/mass) |
| Suitability for iron determination | to pass test |

JIS K 8486:2007

(*:Additional test performed by Wako but not required by JIS)