

Code No.04-0656

SPECIFICATION

3-(3,4-Dihydroxyphenyl)-L-alanine

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown, crystalline powder ~ powder to pass test
Solubility in dil. hydrochloric acid	
Specific rotation [α]D20(c=4,dil.HCl) (after drying)	-10.0~-15.0°
Loss on drying at 105°C	under investigation
Thin layer chromatography test	to pass test
Assay (after drying)	min.98.0%