

Code No.04-0718

## SPECIFICATION

## Dimethyl Adipate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	1.060~1.068g/mL
Freezing point	9~11°C
Refractive index 20°C	1.425~1.432
Water	max.0.2%
Acidity [as HOOC(CH <sub>2</sub> ) <sub>4</sub> COOH]	max.0.1%
Assay (capillary GC)	min.97.0%