

Code No.04-1864

SPECIFICATION

Bis(2-ethylhexyl) Sebacate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.912~0.919g/ml
Refractive index 20°C	1.446~1.454
Water	max.0.1%
Acidity [as HOOC(CH ₂) ₈ COOH]	max.0.02%
Assay (capillary GC)	min.98.0%
Assay (volumetric analysis)	min.95.0%