

SPECIFICATION

N-[3-(Dimethylamino)propyl]acrylamide

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Yellow ~ brownish yellow, clear liquid
Solubility in water	to pass test
Solubility in ethanol	to pass test
Density (20°C)	0.947~0.957g/ml
Refractive index 20°C	1.479~1.485
Water	max.0.5%
Assay (capillary GC)	min.97.0%
Stabilizer (p-methoxyphenol)	0.1~0.5%