

Code No.04-2322

SPECIFICATION

Distilled Water for Preparative Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.997~0.999g/ml
Refractive index 20°C	1.332~1.334
Absorbance (210nm)	max.0.03
Absorbance (220~400nm)	max.0.01
Residue after evaporation	max.5ppm
pH at 25°C	5.0~7.5
Peroxides (as H ₂ O ₂)	max.0.5ppm