

Code No.04-2591

SPECIFICATION

Diethylene Glycol

Wako Special Grade

| REQUIREMENT | SPECIFICATION |
|--------------------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Solubility in water | to pass test |
| Density (20°C) | 1.116~1.123g/ml |
| Water | max.0.1% |
| Acidity (as CH ₃ COOH) | max.0.01% |
| Substances darkened by sulfuric acid | to pass test |
| Assay (capillary GC) | min.99.0% |