

Code No.04-2829

SPECIFICATION

1,2-Dichloroethane
for Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20℃)	1.251~1.258g/mL
Refractive index 20℃	1.442~1.446
Absorbance (235nm)	max.0.15
Absorbance (240nm)	max.0.10
Absorbance (254nm)	max.0.02
Absorbance (270~400nm)	max.0.01
Water	max.0.03%
Fluorescence test	to pass test
Assay (capillary GC)	min.99.9%