

## SPECIFICATION

Diethyl Ether  
for Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.712~0.714g/ml
Absorbance (220nm)	max.0.40
Absorbance (230nm)	max.0.30
Absorbance (254nm)	max.0.10
Absorbance (260nm)	max.0.05
Absorbance (290~400nm)	max.0.01
Water	max.0.05%
Ethanol (C <sub>2</sub> H <sub>5</sub> OH)	max.0.05%
Fluorescence test	to pass test
Assay (capillary GC)	min.99.5%