

SPECIFICATION

Diethyl Ether for Pesticide Residue And Polychlorinated Biphenyl Analysis
Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance Assay (C ₂ H ₅ OC ₂ H ₅) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis * Suitability for pesticide residue analysis (GC-FPD)	Colorless clear liquid min.99.5%(mass/mass) 0.712~0.714g/ml max.0.1%(mass/mass) max.0.0002%(mass/mass) to pass test to pass test

JIS K 8357:2012(*:Additional test performed by Wako but not required by JIS)