

Code No.04-3337

SPECIFICATION

(S)-2,5-Dihydro-3,6-dimethoxy-2-isopropylpyrazine

for Organic Synthesis

| REQUIREMENT | SPECIFICATION |
|--|----------------------------------|
| Appearance | Colorless ~ yellow, clear liquid |
| Solubility in acetone | to pass test |
| Density (20°C) | under investigation |
| Refractive index 20°C | under investigation |
| Specific rotation [α] _{D20} (c=2,CH ₃ OH) | under investigation |
| Water | under investigation |
| Assay (capillary GC) | min.97.0% |