

Code No.04-3397

SPECIFICATION

Dimethyl Sulfoxide

for Dioxins Analysis

| REQUIREMENT | SPECIFICATION |
|---------------------------------------|------------------------|
| Appearance | Colorless clear liquid |
| Density (20°C) | 1.099~1.103g/mL |
| Freezing point | min.18°C |
| Refractive index 20°C | 1.477~1.480 |
| Water | max.0.2% |
| Suitability for dioxins determination | to pass test |
| Assay (capillary GC) | min.99.0% |