

SPECIFICATION

Ethyl Acetate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Description Assay (CH ₃ COOC ₂ H ₅) (GC) Appearance (color) Density (20°C) Refractive index 20°C Water Residue after evaporation Acidity (as CH ₃ COOH) Ethanol (C ₂ H ₅ OH) (GC) Methanol (CH ₃ OH) (GC) Methyl Acetate (CH ₃ COOCH ₃) (GC) Substances darkened by sulfuric acid	Colorless clear liquid min.99.5%(area/area) max.10 Hazen units 0.898~0.902g/mL 1.370~1.374 max.0.1%(mass/mass) max.0.002%(mass/mass) max.0.005%(mass/mass) max.0.1%(area/area) max.0.1%(area/area) max.0.1%(area/area) to pass test

JIS K 8361:2022(*:Additional test performed by Wako but not required by JIS)