

Code No.05-0086

SPECIFICATION

1,2-Dichloroethane

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (CH ₂ ClCH ₂ Cl) (GC) Density (20°C) Water Residue after evaporation Acidity (as HCl) Chloride (Cl) Iodine-consuming substances Substances darkened by sulfuric acid	Colorless clear liquid min.99.5%(mass/mass) 1.251~1.261g/ml max.0.03%(mass/mass) max.0.002%(mass/mass) max.5ppm(mass/mass) max.5ppm(mass/mass) to pass test to pass test

JIS K 8465:2011

*:Additional test performed by Wako but not required by JIS