

Code No.05-0127

## SPECIFICATION

## 2-Ethyl-1,3-hexanediol (mixture of isomers)

Wako 1st Grade

| REQUIREMENT                               | SPECIFICATION                              |
|---|--|
| Appearance                                | Colorless ~ nearly colorless, clear liquid |
| Density (20°C)                            | 0.937~0.945g/ml                            |
| Refractive index 20°C                     | 1.447~1.453                                |
| Water                                     | max.0.1%                                   |
| Assay (capillary GC) (mixture of isomers) | min.98.0%                                  |