

## SPECIFICATION

## Diethyl Malonate

Wako Special Grade

| REQUIREMENT                           | SPECIFICATION                              |
|---------------------------------------|--|
| Appearance                            | Colorless ~ nearly colorless, clear liquid |
| Density (20°C)                        | 1.053~1.061g/ml                            |
| Water                                 | max.0.1%                                   |
| Acidity (as HOOCCH <sub>2</sub> COOH) | max.0.2%                                   |
| Substances darkened by sulfuric acid  | to pass test                               |
| Assay (capillary GC)                  | min.98.0%                                  |