

# SPECIFICATION

## Ethyl Acetate

for Amino Acid Sequence Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (260nm)	max.0.20
Absorbance (270nm)	max.0.05
Absorbance (280nm)	max.0.02
Absorbance (300~400nm)	max.0.01
Water	max.0.02%
Aldehydes (as CH <sub>3</sub> CHO)	max.1 ppm
Peroxides	to pass test
Assay (capillary GC)	min.99.5%