

Code No.05-0380

## SPECIFICATION

Ethyl Lactate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	1.031~1.038g/ml
Refractive index 20°C	1.411~1.415
Specific rotation [ $\alpha$ ]D20(neat)	-0.3~+0.3°
Water	max.0.5%
Acidity [as CH <sub>3</sub> CH(OH)COOH]	max.0.2%
Assay (capillary GC)	min.95.0%