

Code No.05-0419

SPECIFICATION

Ethyl Diethylphosphonoacetate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Solubility in water	to pass test
Solubility in ethanol	to pass test
Density (20°C)	1.120~1.130g/ml
Refractive index 20°C	1.429~1.433
Water	max.0.2%
Acidity (as CH ₃ COOH)	max.0.1%
Assay (capillary GC)	min.97.0%